

Easy Read Hearing Aid Instructions

Information and advice for patients - Audiology

Your hearing aid will look like



The **left**
hearing aid has
a **blue** marker



The **right**
hearing aid has
a **red** marker


Turning the hearing aid on and off

In the morning
close the battery
drawer to turn on
the hearing aid



At bed time
open the
battery drawer
to turn off the
hearing aid

Using the volume control



← Louder



← Quieter

Keeping the hearing aid safe

When you take out
your hearing
aid put it in
the pouch to
keep it safe




Changing the battery

 Open the battery drawer and take out the old battery	 Put in the new battery in and take off the sticker		
 312	 13	 675	 Recycle your old batteries

Batteries are small but very dangerous if swallowed. Keep them safe and seek emergency help immediately if one is swallowed.

Call us

- To get more batteries
- If you are having problems with the hearing aid



0121 507 5366

Audiology, Lyng Health Centre

0121 507 4875

Hearing Services Centre

www.swbh.nhs.uk/services/audiology | Facebook page SWB Audiology

Photographs taken and supplied with permission from Sandwell and West Birmingham NHS Trust, Audiology Department, 2025.

Do you need this information in another language or format? If you require this leaflet in another language or format such as audio or large print, please contact the department or service.

If you would like to suggest any amendments or improvements to this leaflet please contact SWB Library Services on ext 3587 or email swbh.library@nhs.net.